


EXCHANGEABLE COLLOID PROBE

Spherical colloids are the most suitable probes for local elasticity measurements on complex substrates. While colloidal probes are inherently difficult to produce and handle, FluidFM® technology overcomes these limitations in order to give you unparalleled flexibility for your most demanding research requirements.

FluidFM® GIVES YOU THE EDGE.

Imagine renewing your AFM colloidal probe in-situ without having to completely replace the entire probe. FluidFM® technology makes opting for a completely fresh probe inherently easy. The simple, yet universal approach enabled by FluidFM® technology allows it to reversibly attach micro- and nanospheres to an atomic force cantilever in order to function as a colloidal probe.

Quantify long-term or irreversible interactions by using each colloidal probe only once. Fast, in-situ renewal of your probe is possible with FluidFM® tech-



STRONG STATISTICS. 60 µm polystyrene colloids are used to quickly assess cell adhesion. Courtesy of Dörig P., ETH Zurich

nology – at virtually no cost. Obtain solid statistics in short periods of time by measuring more data points than ever before. The versatility of FluidFM® thereby allows you to use solid, liquid and gaseous colloids as required by your experiment.

THE PROCEDURE IN BRIEF.

The colloids are seized and reversibly attached to the FluidFM® probe by applying an underpressure to the microfluidic channel. Once measurements

with the attached colloid concludes, it can be easily detached from the probe by application of a short overpressure pulse.

PUBLICATIONS

- 2013. P. Dörig, D. Ossola, A. M. Truong, M. Graf, F. Stauffer, J. Vörös & T. Zambelli. Exchangeable colloidal AFM probes for the quantification of irreversible and long-term interactions. *Biophysical Journal*, 105(2), 463–472. doi:10.1016/j.bpj.2013.06.002
- 2015. B. R. Simona, L. Hirt, L. Demkó, T. Zambelli, J. Vörös, M. Ehrbar & V. Milleret. Density gradients at hydrogel interfaces for enhanced cell penetration. *Biomater. Sci.* doi:10.1039/C4BM00416G



PROBES PER CANTILEVER

pN

FORCE RESOLUTION

> 1

PROBE PER MINUTE

sub µm
COLLOIDS

CONTACT US.

We offer complete support for our customers and distributors. Please visit the Cytosurge Help Center in order to access the FluidFM® user community. www.fluidfm.com

CYTOSURGE®

CYTOSURGE AG, SÄGEREISTRASSE 25, 8152 GLATTBRUGG, SWITZERLAND
PHONE +41 44 533 14 50, FAX +41 44 533 14 59, WWW.CYTOSURGE.COM

