



Vapor Annealing Effect on Copolymers Studied by Agilent 7500 Atomic Force Microscopy with Environmental Control

Application Brief

Introduction

Microphase separation is a ubiquitous phenomenon in copolymer thin films due to the immiscibility and dissimilar nature of their block components. Such segregations could occur in directions both perpendicular and parallel to the underlying support substrate. In particular, the latter case will expose all of polymer components to the air/polymer interface and form periodic microdomains consisting with distinctive chemical identities. The spatial distribution (i.e.; their geometry, spacing as well as patterning) of the components play an extremely important role in the chemical and physical properties of polymers.

Polymer annealing offers an effective means to tune and regulate those component domains. Generally, annealing can be introduced by heating the copolymer films well above their transition temperatures (T_g) to obtain enhanced microphase separations. However, such thermal annealing could induce some heating-related issues (unexpected polymer oxidation, etc.) and it fuels the research on alternative strategies such as the solvent vapor-assisted annealing. Agilent 7500 AFM system is comprised of a sealed environmental chamber. The desired vapor can be realized readily by putting a small vial of corresponding solution into the chamber. In this application brief,

results on an atomic force microscopy (AFM)-based investigations on the annealing effect on triblock copolymer PS-*b*-PB-*b*-PS (SBS) will be presented.

AFM Structural Characterization of the Triblock Copolymer SBS without Annealing

The initial triblock copolymer SBS samples were prepared by spin-casting of their 10 mg/ml solution in toluene onto a smooth silicon substrate at room temperature and under ambient conditions. Figure 1 are both AFM topography and corresponding phase images of the normal SBS copolymer without any annealing. Compositional mapping with AFM is often used for observations of microphase separation of block copolymers, which occurs at the sub-100 nm scale. This is seen clearly in height and phase images of a triblock copolymer polystyrene-*b*-polybutadiene-*b*-polystyrene (PS-*b*-PB-*b*-PS) film. A microphase separation pattern of this material, which is mostly pronounced in the phase image, is characterized by structural parameter ~ 35 nm. The phase contrast is related to the fact that at room temperature, PS is in a glassy state while PB is in a rubber-like state. Consequently, the brighter areas in the phase image (corresponding higher areas in topographic image) can be attributed to stiff lamellae of PS.




Figure 1. AFM topography and corresponding phase images of copolymer PS-*b*-PB-*b*-PS without annealing. Scan size: $2\mu\text{m} \times 2\mu\text{m}$.



Agilent Technologies





Figure 2. AFM studies of the effect of toluene vapor on the structural changes of SBS copolymer. Scan size: $2\mu\text{m} \times 2\mu\text{m}$.

AFM Structural Characterization of the Triblock Copolymer SBS with Vapor Annealing

When block copolymer films were prepared by the spin-casting method, the resulting product and its surface morphology typically have not reached their equilibrium due to the very fast evaporation of the solvent. This is evident by the fact that those PS domains in Figure 1 are short and discontinuous. The introduction of common solvent vapor such as toluene will lead to a swelling of block copolymers. As a consequence, it could greatly facilitate the molecular motions of components and accelerate the structural evolution towards thermodynamically favored phase segregation. Such process can be monitored in real time by an AFM incorporated with a well-controlled environment. The effect the toluene

vapor on the SBS structures is clearly captured by AFM. Figure 2 (left) exhibits the appearance of long-range ordering of the SBS films after a few hours. As can be seen, the in-plane cylinders become extended and formed the ordered domains as large as hundreds of nanometers. Continuous exposure to the toluene vapor allows the further changes and highly ordered domains at micrometer-scale are observed (Figure 2, right).

Summary

Using investigations of the toluene vapor influence on copolymer SBS surface structures as an example, it is demonstrated that AFM measurements under well-controlled environments can be achieved easily with Agilent 7500 system.

AFM Instrumentation from Agilent Technologies

Agilent Technologies offers high-precision, modular AFM solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Agilent's leading-edge R&D laboratories are dedicated to the timely introduction and optimization of innovative and easy-to-use AFM technologies.

www.agilent.com/find/afm

Americas

Canada	(877) 894 4414
Latin America	305 269 7500
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Thailand	1 800 226 008

Europe & Middle East

Austria	43 (0) 1 360 277 1571
Belgium	32 (0) 2 404 93 40
Denmark	45 70 13 15 15
Finland	358 (0) 10 855 2100
France	0825 010 700*
	*0.125 €/minute
Germany	49 (0) 7031 464 6333
Ireland	1890 924 204
Israel	972-3-9288-504/544
Italy	39 02 92 60 8484
Netherlands	31 (0) 20 547 2111
Spain	34 (91) 631 3300
Sweden	0200-88 22 55
Switzerland	0800 80 53 53
United Kingdom	44 (0) 118 9276201

Other European Countries:

www.agilent.com/find/contactus

Product specifications and descriptions in this document subject to change without notice.

© Agilent Technologies, Inc. 2013
Published in USA, September 16, 2013
5991-3251EN



Agilent Technologies