

# H4KFWSD08DPX

- ▶ 4K resolution
- ▶ On screen measurement
- ▶ 4K Video record
- ▶ Digital Classroom with APP
- ► HDMI, USB, GE and WiFi



















### H4KFWSD08DPX

The H4KFWSD08DPX produces high-resolution images with vibrant natural colors suited also for critical applications like Pathology and likewise plus it has a state of the art built-in measurement software providing highly effective measurement tools directly on the screen with no need for a computer.

With the H4KFWSD08DPX users can enjoy fast, comfortable, and user-friendly on-screen operation — improving sample evaluation, measurement, facilitating effective collaborations, and documentation in high resolution. For more advanced applications the user has the option to connect the camera to a PC via USB 3.0, GE, or WiFi and use DeltaPix InSight in order to make further analysis. Moreover, this camera provides a completely new platform for educational applications because it can stream the video to android smartphones making it super easy to share images and videos with the students



#### High quality, cost effective and flexible

The DeltaPix H4KFWSD08DPX microscope camera has C-mount interface making them easily and cost-efficiently fitted on most standard microscopes. Additionally, it has multiple interfaces providing unmatched versatility for any use case.



## User-friendly Camera

Plug and play — the HDMI interface provides high speed and ease of installation on any 4K montitor with HDMI input.

No need for computer, because the H4KWSD08DPX microscope camera has a built-in software. Option to use it with PC software of Smartphone app.




**Typical applications:** ✓ Documentation and publishing ✓ Material science ✓ Metrology ✓ Quality control ✓ Dark Field ✓ Bright Field ✓ Histology ✓ Pathology ✓ Semiconductor inspection ✓ Digital Classroom

Spectral Sensitiviy

#### Spectral Sensitivity Characteristics


(Excludes lens characteristics and light source characteristics.)



# Why choose H4KFWSD08DPX?

#### DeltaPix Build-in Software

The H4KFWSD08DPX camera has a the stat of art software built-in software that provides a powerful platform for precise measurements, image capturing, and video recording.



#### **Onboard measurement tools:**

- Angle
- Four point angle
- Point
- Arbitrary line
- Three point line
- Horizontal line
- Vertical line
- Parallel
- Three point vertical
- Rectangle
- Curve

- Ellipse
- Five point ellipse
- Three pint circle
- Center + radius circle
- Arc
- Text
- Polygon
- Annulus





#### **Camera Dimensions**



# Features

#### H4KFWSD08DPX H4KFWSD08DPX 3840 x 2160 @60FPS (HDMI) Live preview 1920 x 1080 @60FPS (HDMI) 3840 x 2160 @30FPS (USB 3.0) 3840 x 2160 @30FPS (WiFi) 8.3Mpixels (3840x2160) Record video 4K Video Built-in, UI in English /Chinese / Korean / Thai / French / German / Japanese / Italian InSight Basic (Included) \*Free Android APP 1/1.8" Sony back illuminated Exmor sensors 7.68mm x 4.32mm 2μm x 2μm G Sensitivity 505mv with 1/30s 0.1mv with 1/30s Automatic or Manual Automatic, Manual or ROI HDMI2.0, USB 3.0, WiFi and GE Still image, JPEG or Tiff H264 encoded video Minimum 4K/Full HD monitor with HDMI input

Mouse, 32GB High speed SD card and 12V power supply.

#### \*Coming soon

DeltaPix Aps, Hassellunden 16, DK-2765 Smorum Denmark Telephone:+45 46760205 E-mail: deltapix@deltapix.dk